




# **MEDPOR® Biomaterials CRANIAL/NEUROSURGICAL**

Innovative Technology in Reconstructive Surgical Implants



# MEDPOR® Customized and Off-the-Shelf-Implants — Solutions for Neurosurgeons



MEDPOR® Porous Polyethylene Implants provide surgeons with an expanding range of options for reconstruction and augmentation. Now, Porex Surgical offers surgeons more options to correct cranial defects:

**Standard Customized Implant Services** – For surgeons who prefer a physical skull model and template, or for some complex bilateral defects, or defects involving the orbital or facial structures.

**Customized Implant Services with e-viewCT™** – The fastest, easiest way to order, view and obtain a MEDPOR Customized Implant is via a virtual process via Porex Surgical's secure Web utility.\*

**Off-the-Shelf Implants for Complex Procedures** – Porex Surgical's off-the-shelf implants for cranial, neurosurgical and skull-base applications provide surgeons with an attractive alternative to bone graft procedures.

\*Due to the size and complexity of some defects, some customized implants may need to be produced in more than one piece, may not be able to be viewed on-line, and may require a physical model to be produced and shipped for review.

## Surgical Planning is Faster and Easier with MEDPOR® Customized Implants

MEDPOR Customized Implants provide the surgeon with an attractive alternative to complex grafts and other implant materials. Porex Surgical Customized Implant Services can provide implants shaped for a patient's individual needs. A 3-D model of the patient's bony anatomy is created from CT scan data; Customized Implant Shapes can then be made to fit the defect or to correct an asymmetry. MEDPOR Customized Implant Services are available as standard services or with **e-viewCT™**.

**Standard Customized Implant Services** – MEDPOR Customized Implants for complex bilateral defects, or defects involving the orbital or facial structures are created using a physical skull model and template.



**Customized Implant Services with e-viewCT** – The fastest, easiest way to order, view and obtain a MEDPOR Customized Implant is via Porex Surgical's new and improved customized implant secure Web application.

**Advantages:**

- Ability to transmit the patient's CT scan data electronically via a secure Web portal
- High quality, customized implants that may significantly reduce operating time and expense


**Steps to Creating MEDPOR Customized Implants**

Creating a MEDPOR Customized Implant for an individual patient's complex bilateral defect or defects involving the orbital floor is a multi-step process and requires close communication between the surgeon and Porex Surgical.

- CT data should be obtained using Porex Surgical's Scanning Protocol.
- A purchase order is submitted to Porex Surgical. The customized implant process cannot be started without receipt of a purchase order or pre-payment of the implant.
- CT data is submitted to Porex Surgical via the secure **e-viewCT** utility or on a CD.
- Porex Surgical converts the patients CT scan data into either a virtual model and customized implant template that can be viewed and rotated 360 degrees electronically via **e-viewCT** or a physical model and non-sterile implant template.
- A prescription form is provided to the surgeon for review and approval of the final implant shape.
- Upon receipt of the signed prescription form signifying physician approval of the template shape, a MEDPOR Customized Implant is manufactured, sterilized, and shipped.



***Feathering of the Implant***



***Fixation of Implant***



***Immediate Pre-Op***



***Twenty-One Months Post-Op***

Photos courtesy of Robert D. Wallace, MD

# e-viewCT™

advanced technology for innovative surgeons



## New and Improved!

Porex Surgical's new and improved interactive Web portal allows for faster and easier transfer of patient data between our highly skilled design staff and the surgeon or radiologist. Customized implant shapes can then be made to fit the defect or to correct an asymmetry.\* From receipt of CT scan data, to a virtual rotating 3-D skull model and implant template design, to the shipment of the final sterile implant, the redesigned **e-viewCT** Web tool allows the surgeon 24-hour access to the process. The improved tool allows the surgeon to log in with secure access and view the status of a specific case. An individual date stamp is documented as each of the milestones that make up the customized implant process is completed. E-mail notifications of the process can be automatically generated to individuals designated by the surgeon, enabling more efficient planning through faster communication. **e-viewCT** is advanced technology for innovative surgeons.

\*Due to the size and complexity of some defects, it may be necessary to produce some customized implants in more than one piece. It may not be possible to view these implants on line and it may be necessary to produce a physical model and ship it to the surgeon for review.

## MEDPOR® Customized Cranial Implant



CAT#	DESCRIPTION
89020*	MEDPOR Customized Cranial Implant (Includes Peri-Orbital Customized Implants) Package Includes: <ul style="list-style-type: none"> <li>• On-line review of skull model (if necessary)</li> <li>• On-line approval of implant template (if necessary)</li> <li>• One sterile customized implant plus one sterile backup customized implant (USA only)</li> </ul>
89024**	Cranial implant template add-on, non-implantable

Delivery time is approximated from receipt of purchase order and CT data at Porex Surgical. Call for an estimated delivery time. Note: Complex bilateral defects, or defects involving the orbital structures or facial structures, may require additional expense, time, and/or a physical skull model and template. Please call for a quote.

## MEDPOR® Customized Facial Implants



CAT#	DESCRIPTION
89021*	MEDPOR Customized Facial Implant (Includes Chin, Mandible, Malar, and Midface Implants) Package Includes: <ul style="list-style-type: none"> <li>• One sterile customized implant plus one sterile backup implant (USA only)</li> </ul>
89022	Contralateral charge for Customized Facial Implant
89023**	Skeletal model of defect area add-on
89025**	Facial implant template add-on, non-implantable
89026**	Bilateral implant template add-on, non-implantable

\*Catalog #'s 89020 & 89021 DO NOT include a physical model or template.




\*\*Model price is valid only with purchase of a MEDPOR Customized Implant



# MEDPOR Biomaterial

Both MEDPOR® Customized Implants and off-the-shelf implants are manufactured from MEDPOR Biomaterial. MEDPOR Porous Polyethylene Implants provide surgeons with an expanding range of options for reconstruction and augmentation. MEDPOR is a biocompatible, porous polyethylene material. The interconnecting, omni-directional pore structure allows for fibrovascular in-growth and integration of the patient's tissue. More than 250,000 procedures have been performed with MEDPOR Biomaterial, with more than 350 published clinical reports in cranial, reconstructive, oculoplastic and cosmetic applications.


- MEDPOR is easy to work with! The material can be trimmed with a blade in the sterile field, carved and feathered intra-operatively for an excellent final fit.
- No pre-placing of fixation plates! MEDPOR can be easily drilled and fixated and will accept screws and plates without cracking, giving the surgeon more flexibility in fixation options and placement.
- MEDPOR Surgical Implants are easily cut with a variety of surgical instruments. The physical properties of the MEDPOR Biomaterial allow for cutting and trimming the implant while maintaining the interconnectivity and the structure of the pores. Implants may require fitting to the defect area at the time of surgery. The implant edges should be delicately shaped and feathered for a smooth transition from the implant to the patient's own bony contour.
- MEDPOR Surgical Implants are provided sterile and should not be resterilized.
- Do not place or carve the implant on surgical drapes, surgical clothing or any other surface that may contaminate the implant with lint and other particulate matter.



1. Sheets and blocks trim easily with surgical scissors or a small scalpel blade.
2. Thicker implants may be shaped with large scalpels, bone cutters, or a cutting burr.
3. Multiple pieces of implant material may be stacked and sutured together.
4. Proper rigid fixation techniques allow for stabilization of the implant.
5. Submerge the implant in a hot, sterile saline bath (over 180°F) for several minutes until the implant's memory relaxes.
6. Gently bend the implant to the desired shape. Return the implant to the hot saline if there is too much resistance.
7. Hold the implant in the final shape and allow to cool. The use of a cold sterile bath can accelerate this process.\*

\*Thicker or larger MEDPOR Implants may be difficult to shape using this method.

## MEDPOR TITAN® IMPLANTS



### Strength Meets Flexibility

*Designed with: Nicholas T. Iliff, M.D., Shannath L. Merbs, M.D., Ph. D., and Michael P. Grant, M.D., Ph. D.*

MEDPOR TITAN® Sheets are intended for non-weight bearing applications of craniofacial reconstructive/cosmetic surgery, and repair of craniofacial trauma.

Titanium mesh and MEDPOR Polyethylene Implants have a long history of successful use in trauma repair. When cut, traditional titanium mesh may exhibit many sharp points and edges that can make insertion difficult. Titanium mesh embedded within a thin sheet of high-density polyethylene may minimize sharp edges even when the implant is cut. The titanium mesh is radiopaque, making the implant visible on radiographs or CT scans. The titanium mesh used in MEDPOR Biomaterial allows the surgeon to bend and contour a thin implant material to the desired shape while providing the strength usually associated with a much thicker traditional MEDPOR Implant.

U.S. Patent #7,655,047


中国发明专利, 专利号ZL200480009959.6

Surgeons may choose from three types of MEDPOR TITAN Sheets:

- The MEDPOR TITAN MEDPOR (MTM™) Implant is porous, high-density polyethylene with titanium mesh embedded in it, providing the advantages of fibrovascular integration of the patient's host tissue through the sheet.
- The MEDPOR TITAN BARRIER (MTB™) Implant is a sheet of titanium mesh embedded within a porous polyethylene matrix with a solid, BARRIER™ surface on one side, allowing for fibrovascular ingrowth only on the porous side of the implant.
- The MEDPOR TITAN Double BARRIER (BTB™) Implant is titanium mesh embedded within solid, high-density polyethylene that acts as a BARRIER to tissue attachment and may help facilitate implant placement.

CAT#	DESCRIPTION	A	B	C
81020	MTM	76mm x	50mm x	0.85mm
81021	MTM	38mm x	50mm x	0.85mm
81022	MTM	38mm x	50mm x	1.5mm
81023	MTM	76mm x	50mm x	1.5mm
81024	BTB	38mm x	50mm x	0.6mm
81025	BTB	76mm x	50mm x	0.6mm
81026	MTB	38mm x	50mm x	1.0mm
81027	MTB	76mm x	50mm x	1.0mm
81028	MTB	38mm x	50mm x	1.6mm
81029	MTB	76mm x	50mm x	1.6mm

## MEDPOR TITAN® CRANIAL CURVE




The MEDPOR TITAN® Cranial Curve and MEDPOR TITAN® Cranial Curve - BARRIER™ Implants offer the cranial surgeon an attractive option for cranial/skull base reconstruction. Both configurations are pre-shaped to the general curvature of the cranium and are intended for non-loading bearing applications.

The titanium mesh used with the MEDPOR Biomaterial allows the surgeon to further bend and contour the implant material to the desired shape.

Titanium is radiopaque, making the implant visible on radiographs or CT scans.

U.S. Patent #7,655,047


中国发明专利, 专利号ZL200480009959.6



### MEDPOR TITAN Cranial Curve Implant

CAT#	DESCRIPTION	A	B	C
82019	TITAN Cranial Curve	38mm x	48mm x	0.85mm
82020	TITAN Cranial Curve - BARRIER	38mm x	48mm x	1.00mm

## MEDPOR TITAN® POSTERIOR IMPLANT



**NEW !**


*Designed with: Jason Sheehan, MD, PhD, FACS*

The MEDPOR TITAN® Posterior Implant is intended for reconstruction of the cranium.

- Provides an option to other repair materials for reconstructing the cranium.
- May be trimmed and cut with surgical scissors.
- Polyethylene coating may minimize sharp edges of titanium when cut.
- Titanium mesh embedded in the MEDPOR Biomaterial helps the implant retain its shape when bent and contoured to meet a specific patient defect.

U.S. Patent #7,655,047

中国发明专利, 专利号ZL200480009959.6



### MEDPOR TITAN Posterior Implant

CAT #	DESCRIPTION	A	B	C
82030	Posterior Implant	96mm x	61mm x	1.5mm

MEDPOR TITAN® MAX™ SHEET



The new MEDPOR TITAN® MAX™ Sheet is intended for non-weight-bearing applications of craniofacial reconstruction and repair of craniofacial trauma where a larger length and width implant is desired.

The MEDPOR TITAN MAX Sheet is an excellent option to bare titanium mesh for general cranial repair of small to medium-sized defects. The titanium mesh used in the MEDPOR Biomaterial helps the implant retain its shape, which allows the surgeon to bend and contour the implant material to fit a patient-specific defect.


Provides MAX Options for Craniofacial Reconstruction

- Wider titanium mesh in a thin sheet — 1.5mm sheet thickness
- Larger length and width — 76mm x 100mm
- Excellent option to bare titanium mesh for general cranial repair of small to medium size defects

Provides MAX Benefits

- Biocompatible MEDPOR — 20 years of proven use in CMF applications — low risk for complications
- Easily shaped and cut — convenient, may save time, fit to individual patient contours
- Easily fixated with plates/screws - stays in place
- Allows for tissue ingrowth — enhances stabilization and possible reduced risk of long-term complications
- Polyethylene coating may minimize sharp edges of titanium when cut
- Titanium mesh is radiopaque — visible on postoperative radiographs and CT scans

U.S. Patent #7,655,047  
中国发明专利， 专利号ZL200480009959.6



MEDPOR TITAN MAX

CAT #	DESCRIPTION	A	B	C
81040	MTM	100mm x	76mm x	1.5mm


MEDPOR TITAN® CRANIAL-TEMPORAL



The MEDPOR TITAN Cranial — Temporal Implant is designed to be an off-the shelf solution for non-weight-bearing applications of craniofacial reconstruction and repair of craniofacial trauma. The titanium mesh embedded in the MEDPOR Biomaterial provides strength usually associated with a much thicker traditional MEDPOR Implant. The thinner profile can be bent to the shape of the defect. The radio-opaque titanium mesh makes the implant visible on postoperative radiographs or CT scans and helps the implant retain the shape when bent. The implant is available in left and right configurations.

Sterile template included.

U.S. Patent #7,655,047  
中国发明专利， 专利号ZL200480009959.6



MEDPOR TITAN Cranial-Temporal

CAT#	DESCRIPTION	A	B	C
81037	TITAN Cranial - Temporal, Left	130mm x	130mm x	3mm
81038	TITAN Cranial - Temporal, Right	130mm x	130mm x	3mm




## MEDPOR® CRANIAL HEMISPHERE



The MEDPOR® Cranial Hemisphere for large cranial defects provides surgeons with an off-the-shelf alternative to customized implants, complex grafts, and other implant materials.

The implant shape approximates the contour of the half cranium. The Cranial Hemisphere is available in two thicknesses and left and right versions. The implant can be trimmed to fit the defect with a blade. The edges should be delicately shaped and feathered with surgical scissors or a scalpel blade for a smooth transition from the implant to the patient's bony contour. Fixation of the implant may be accomplished with suture, surgical wire, or craniofacial rigid fixation plates and screws.




### MEDPOR Cranial Hemisphere

CAT#	DESCRIPTION	A	B	C
82000	Cranial Hemisphere - Right	124mm	170mm	4.5mm
82001	Cranial Hemisphere - Left	124mm	170mm	4.5mm
82002	Cranial Hemisphere - Right	133mm	170mm	6mm
82003	Cranial Hemisphere - Left	133mm	170mm	6mm

## MEDPOR® CRANIAL DOME



The MEDPOR® Cranial Dome is designed to provide surgeons with a reconstructive option for large cranial defects that encompass the dome area of the skull. The Cranial Dome approximates the contour of the superior one-third of the cranium and is available in two thicknesses.



### MEDPOR Cranial Dome


CAT#	DESCRIPTION	A	B	C
82004	Cranial Dome	150mm	180mm	4mm
82006	Cranial Dome	150mm	180mm	6mm

## FLEXBLOCK™ CRANIAL GRID



The FLEXBLOCK Cranial Grid is designed to fill full thickness cranial defects as an option to calvarial bone grafts. The Cranial Grid has a grid design on the interior surface that provides strength and flexibility and allows the implant to be cut to the desired shape. The shape approximates the contour of the cranium and may be customized to the desired shape by relaxing the memory and bending the implant after submersion in hot, sterile saline. The implant is 6mm thick and 97mm x 106mm in dimension.

U.S. Patent #5,545,226



### MEDPOR FLEXBLOCK Cranial Grid Implant

CAT#	DESCRIPTION	A	B	C
9524	Cranial Grid	97mm	106mm	6mm




FLEXBLOCK™ IMPLANT



The FLEXBLOCK™ Implant is available for small to medium sized cranial defects and contour deformities. It can be used as an onlay for small calvarial defects and contour defects.

The FLEXBLOCK Implant has a smooth exterior surface and a series of pedicles on the interior surface that provide volume and flexibility. The outer perimeter and corresponding pedicles of the implant should be trimmed to provide a flange just larger than the defect for support and fixation.

U.S. Patent #5,545,226



- 1. The pedicles on the inferior surface of the implant provide flexibility to fit the implant to the desired contour.
- 2. The implant should be trimmed just larger than the defect. After cutting the implant to shape, the pedicles along the perimeter should be trimmed to provide a flange for fixation to the underlying bone.
- 3. A “shelf” in the surrounding bone may be created at the edge of the defect to provide for a more smooth transition.

MEDPOR FLEXBLOCK Implant

CAT#	DESCRIPTION	A	B	C
6314	FLEXBLOCK	56mm x 91mm	x 4mm	
82022	FLEXBLOCK - BARRIER	56mm x 91mm	x 4mm	


FLEXBLOCK™ TF & TF2



The FLEXBLOCK TF2 and FLEXBLOCK TF Implant shapes are designed to augment deficient soft tissue in the temporal region.

The FLEXBLOCK TF2 comes in left and right versions and has a thinner, contoured temporal surface for an improved anatomical fit and a more natural post-op result.

U.S. Patent #5,545,226



MEDPOR FLEXBLOCK TF2 & TF

CAT#	DESCRIPTION	A	B	C
9857	TF2, Small - Left	61mm x 78mm	x 18mm	
9858	TF2, Small - Right	61mm x 78mm	x 18mm	
9859	TF2, Med. - Left	74mm x 93mm	x 20mm	
9860	TF2, Med. - Right	74mm x 93mm	x 20mm	
9861	TF2, Large - Left	82mm x 105mm	x 20mm	
9862	TF2, Large - Right	82mm x 105mm	x 20mm	
9521	TF Small	70mm x 70mm	x 10mm	
9522	TF Med.	88mm x 86mm	x 15mm	
9523	TF Large	98mm x 95mm	x 18mm	

PTERIONAL IMPLANT




Designed with: Murali Guthikonda, M.D., F.A.C.S.

The Pterional Implant is designed to correct temporal hollowing in patients who have had surgery involving the pterional approach to the brain. While the pterional craniotomy is one of the most versatile approaches in neurosurgery, it can lead to temporal hollowing. The implant is placed deep to the temporalis during closure to correct this defect.

The Pterional Implant is available in left and right versions, and although similar to the FLEXBLOCK™ TF, is much smaller in design to provide appropriate augmentation.


U.S. Patent #5,545,226



MEDPOR PTERIONAL Implant

CAT#	DESCRIPTION	A	B	C
9864	Right	44mm x 43mm	x 6mm	
9865	Left	44mm x 43mm	x 6mm	

## MEDPOR SHEETS AND BLOCKS



MEDPOR Biomaterial Sheets provide the surgeon with excellent options for craniofacial reconstruction and augmentation. The individually packaged, sterile implants provide “off-the-shelf” availability, and may save time and the expense of harvesting graft material. MEDPOR Sheets are available in a variety of sizes and in thicknesses ranging from 0.25mm to 3.0mm. Feathering the edge of the sheets allows for a smooth transition from the implant to the adjoining skeletal structure.

### MEDPOR Micro Thin Sheets


CAT #	A		B		C
83020	38mm	x	50mm	x	0.25mm
83021	76mm	x	50mm	x	0.25mm
83022	38mm	x	50mm	x	0.35mm
83023	76mm	x	50mm	x	0.35mm
8438	30mm	x	50mm	x	0.40mm
83029	38mm	x	50mm	x	0.45mm
83030	76mm	x	50mm	x	0.45mm

### MEDPOR Ultra Thin Sheets

CAT #	A		B		C
7210	38mm	x	50mm	x	0.85mm
7212	50mm	x	76mm	x	0.85mm
7214	76mm	x	127mm	x	0.85mm
7216	127mm	x	178mm	x	0.85mm

### MEDPOR Sheets

CAT #	A		B		C
6330	38mm	x	50mm	x	1.5mm
6331	50mm	x	76mm	x	1.5mm
8662	76mm	x	127mm	x	1.5mm
6351	127mm	x	178mm	x	1.5mm
9562	38mm	x	50mm	x	3.0mm



The surgeon can carve thicker implants in the sterile O.R. field to obtain implant contours individualized for the surgical situation. Allowing the implant to soak several minutes in a hot, sterile saline bath will relax the memory of the implant, enabling modification of the shape.

### MEDPOR Blocks

CAT #	A		B		C
6332	13mm	x	38mm	x	3mm
6333	13mm	x	38mm	x	6mm
6334	13mm	x	38mm	x	9.5mm
6335	25mm	x	50mm	x	3mm
6336	25mm	x	50mm	x	6mm
6337	25mm	x	50mm	x	9.5mm
6338	38mm	x	63mm	x	3mm
6339	38mm	x	63mm	x	6mm
6340	38mm	x	63mm	x	9.5mm

## MEDPOR TSI™




*Designed with: Murali Guthikonda, M.D., F.A.C.S.*

The MEDPOR TSI™ is designed to repair the sellar floor.

The TSI Implant is available in two sizes and configurations. The larger TSI is designed with three small tabs oriented to facilitate handling and placement while the original TSI design has a single tab.

A nonporous sheet of polyethylene heat-bonded to the posterior surface of the larger TSI forms a BARRIER to aid in preventing tissue ingrowth.



### MEDPOR TSI™

CAT#	DESCRIPTION	A	B	C	D
82007	TSI	20mm	x 20mm	x 2.5mm	x 0.45mm
82008	TSI BARRIER - Large	40mm	x 40mm	x 2.5mm	x 0.73mm



## MASTOID IMPLANT




*Designed with: Dennis Bojrab, M.D.*

For patients undergoing cranial procedures that require removal of bone in the mastoid area, the MEDPOR Mastoid Implant provides surgeons with a convenient method to repair defect areas.

The implants are available in small and regular sizes and should be trimmed at the time of surgery to fit the needs of the individual patient, allowing edges to overlap the defect area by approximately 0.5cm.

The regular Mastoid Implant is available in left and right configuration, while the small Mastoid Implant provides a universal fit to either the left or right side.



### MEDPOR Mastoid Implant


CAT#	DESCRIPTION	A	B	C
82014	Mastoid Small	36mm	x 45mm	x 1.0mm
82015	Mastoid Implant - Left	58mm	x 56mm	x 1.5mm
82016	Mastoid Implant - Right	58mm	x 56mm	x 1.5mm

## ORBITO-ZYGOMATIC (OZ™)



*Designed with: Saleem Abdulrauf, M.D.*

The MEDPOR Orbito-Zygomatic (OZ™) Implant is designed for reconstruction of the superior and lateral surfaces of the orbital roof. The MEDPOR OZ Implant provides surgeons with a convenient "off-the-shelf" anatomically shaped implant to cover the bony orbital roof and lateral wall removed during cranial procedures. The OZ Implant, available in left and right versions, should be trimmed at the time of surgery to fit the needs of the individual patient.




### MEDPOR Orbito-Zygomatic (OZ)

CAT #	DESCRIPTION	A	B	C
81013	Left	33mm	x 38mm	x 0.8mm
81014	Right	33mm	x 38mm	x 0.8mm


## COMPLETE & 2/3 ORBIT SHAPES



Complete and Inferior 2/3 Orbit Implants are designed to replace non-load bearing, bony structures of the orbital area. Complete and 2/3 Orbits are typically carved with a blade, scissors or burr to fit the patient's defect and fixed with sutures, wires or craniofacial screws and plates.



**Complete Orbit**



**Inferior 2/3 Orbit**

### MEDPOR Complete & 2/3 Orbit Shapes


CAT #	DESCRIPTION	A	B
9567	Inferior 2/3 Orbit - Left	108mm	x 75mm
9568	Inferior 2/3 Orbit - Right	108mm	x 75mm
9569	Complete Orbit - Left	93mm	x 75mm
9570	Complete Orbit - Right	93mm	x 75mm

## CRANIOTOMY GAP WEDGE V & T



*Designed with: Professor Jaechan Park, M.D.*

The MEDPOR Craniotomy Gap Wedge V & T are designed to fill a gap often left after a craniotomy. The Craniotomy Gap Wedge V & T implants are triangle shaped implants designed to fit snugly into the gap along a bone flap. The revised Craniotomy Gap Wedge V Implant is a wedge-shaped strip that has a "U" shaped cross section. The Craniotomy Gap Wedge T Implant is designed with a thin flat section on the top surface extending 3mm on each side, for a total roof width of 10mm. Both designs measure 102mm in length. The implants are sold one of each style for a total of two implants per sterile package.



### MEDPOR Craniotomy Gap Wedge V & T


CAT#	DESCRIPTION	A	B	C
82011	Craniotomy Gap Wedge V (1 per package)	102mm	4mm	3.6mm
	Craniotomy Gap Wedge T (1 per package)	102mm	6mm	10mm

## CRANIOTOMY GAP FILLER



*Designed with: Murali Gurthikonda, M.D., F.A.C.S.*

The MEDPOR Craniotomy Gap Filler is designed to fill a gap often left after a craniotomy. The Craniotomy Gap Filler has a flat top surface to cover over the width of the gap with a thin "V" shaped extension designed to be pushed into the gap and expand.




### MEDPOR Craniotomy Gap Filler

CAT#	DESCRIPTION	A	B	C	D
82013	Craniotomy Gap Filler (4 per package)	30mm	6mm	10mm	0.6mm

## BURR HOLE COVERS



Burr Hole Covers are designed to fit into and over holes made by a cranial perforator. The large cover has a 14mm diameter lotus type stem allowing easy size modification. The smaller cover is designed for 5mm diameter holes. The superior flange covers any gaps between the cranial hole and the bone flap.



### MEDPOR Burr Hole Covers

CAT#	DESCRIPTION	A	B	C
7511	Burr Hole Cover	15mm	3mm	5mm
7512	Burr Hole Cover	29mm	7mm	14mm

**POREX**  
SURGICAL PRODUCTS GROUP

Please contact Porex Surgical, Inc. for a complete list of craniofacial and reconstructive MEDPOR Shapes. E-mail: [surgical.info@porex.com](mailto:surgical.info@porex.com)

© 2010 Porex Surgical, Inc. MEDPOR and MEDPOR TITAN are registered trademarks and FLEXBLOCK, TSI, e-viewCT, OZ, MTM, BTB and MTB are trademarks of Porex Surgical, Inc. The POREX SURGICAL service mark is owned by or used under the license and authority of Porex Corporation.

#### Porex Surgical, Inc.

15 Dart Road, Newnan, GA 30265-1017 USA  
Email: [surgical.info@porex.com](mailto:surgical.info@porex.com)  
Telephone: 1-678-479-1610 • 1-800-521-7321  
Fax: 1-678-423-1435 • [www.porexurgical.com](http://www.porexurgical.com)

#### Porex Surgical GmbH

Balanstrasse 55, 81541 München, Germany  
Email: [surgical.info2@porexurgical.com](mailto:surgical.info2@porexurgical.com)  
Telephone: +49 (89) 232415-12  
Fax: +49 (89) 232415-15 • [www.porexurgical.de](http://www.porexurgical.de)