SURGICAL PRODUCTS GROUP

**MEDPOR® SST™, SST™ EZ™,
SST™ MCOI™, PLUS™ SST™ and
PLUS™ SST™ EZ™ Orbital Implants**

Suggestions for Use

Porex Surgical Inc.

15 Dart Road Newnan, GA 30265-1017 USA

1-678-479-1610 / 1-800-521-7321 / Fax: 1-678-423-1435

www.porexsurgical.com

MEDPOR® Surgical Implants

Suggestions for Using the MEDPOR® SST™, SST™ EZ™, SST™ MCOI™, PLUS™ SST™ and PLUS™ SST™ EZ™ Orbital Implants

Shape Design:

The MEDPOR® SST™, SST™ EZ™, SST™ MCOI™, PLUS™ SST™ and PLUS™ SST™ EZ™ Orbital Implants have a smooth porous anterior surface and four shallow suture tunnels on the anterior surface that facilitate suturing of the rectus muscles to the implant. These spherical implants are designed to fill the void left by the eviscerated or enucleated globe. They are manufactured from high density polyethylene which has been shown to have good tissue tolerance and stability. The smoother porous surface may help to minimize late implant exposures.

Tunnel Location:

Four tunnels are arranged in a square pattern around the anterior apex of the implant. From an anterior apex view the tunnels connect two adjacent (neighboring) holes (Fig. 1).

The tunnels run from one hole to the next adjacent one and are curved for easy insertion of ophthalmic needles typically used to suture rectus muscles to orbital implants (Fig. 2). The suture needle should follow this same track.

The surgeon should be thoroughly familiar with the information contained in the MEDPOR Surgical Implant package insert prior to use of any MEDPOR Surgical Implant.

Applications for Use:

These shapes are designed for use after standard enucleation or evisceration procedures. The implants are primarily designed (1) to fill the void resulting from an enucleated or eviscerated eye, (2) with a smooth porous anterior surface which may minimize late implant exposures, and (3) to allow easy rectus muscle attachment directly to the MEDPOR Biomaterial.

Contraindications:

Please see the contraindications in the MEDPOR and MEDPOR PLUS Product Information Sheets included in the product package.


Fig. 1


Fig. 2

Shaping the Implant

Good surgical judgment must be used in determining the appropriate size and shape of the final implant for each patient. The porous implant should not be placed on surgical drapes or other materials that may adhere to its textured surface. During surgery, it may be convenient to place the implant on the inner sterile pouch, or in a basin of sterile saline.

If the implant is modified at the time of surgery, sharp or rough edges of the implant should be feathered with a scalpel, and there should be no irregularities present on the implant.

Prior to implantation, the MEDPOR Surgical Implant may be rinsed in a sterile antibiotic solution.

In the case of MEDPOR PLUS SST or MEDPOR PLUS SST-EZ, the implant should be saturated with sterile saline before implantation. Refer to the MEDPOR PLUS Product Information Sheet included in the product package for detailed instructions.

Operative Technique:




In enucleation, the globe is removed in the usual fashion. The implant is oriented within the intraconal compartment, so that the smooth porous surface with the suture tunnels is oriented anteriorly facing away from the posterior apex of the orbit. The tunnels are designed to accept the suture needles commonly used to attach the extra-ocular muscles to the implant.

The needle of the suture should be passed through the tunnel between two adjacent holes (Fig. 3, Fig. 1). Use a fine needle holder and grasp the hub of the needle near its attachment to the suture. Carefully pass the needle through the tunnel following the path of the tunnel as it traverses the implant. The curve of the needle should allow for easy passage through the tunnel. Both needles of the double-armed suture may be passed through the tunnels in a "reverse" or "crossed swords" fashion (Fig. 4).

An alternative suture placement is one in which one arm of the double-armed suture is passed through the tunnel and then tied to the remaining arm of the suture. If a "crossed swords" technique is used, one must be careful when passing the second needle not to cut the suture already inside the tunnel.

The implant is designed to allow the surgeon discretion in muscle placement.

If the muscles are snugged tightly against the implant suture holes by full tightening of the sutures, the muscles will rest at an anterior position on the implant (Fig. 5).

If the surgeon prefers a more equatorial attachment of the muscle against the implant, the sutures can be tied more loosely allowing the muscle to "hang back". By adjusting the tension of the sutures, each of the rectus muscles can be positioned at the chosen point of the implant.

If the surgeon desires, a stronger needle may be used to penetrate the implant, thus attaching the muscles to the implant in any desired position. Placement of a suture into the more dense anterior area of the implant may require a larger or more stout needle than the needle commonly used to attach extra-ocular muscles. Once muscle attachment has been accomplished with all remaining rectus muscles, Tenon's capsule and conjunctiva are closed in separate layers.

In evisceration procedures, posterior relaxing incisions are made in the scleral shell to provide adequate space for the implant. This is particularly important when evisceration is performed in small or phthisical eyes. The implant is then placed within the scleral shell with the smoother surface and suture tunnels directed anteriorly as described for enucleation. The sclera is closed over the implant, and Tenon's capsule and conjunctiva are closed in the usual fashion.

Additional Suture Placement Suggestions:


Suggestion #1

When a suture needle is placed in a needle holder in the normal fashion with the point at a right angle (Fig. 6), following the path of the tunnel in situ may become difficult (Fig. 7), as the needle holder is angled to enter the orbit.

Grasping and rotating the implant (Fig. 8) with toothed forceps so that the tunnel is near the anterior apex of the orbit will tilt the tunnel path so that it approximates the direction of the needle. A curved rather than a straight needle holder may facilitate passage of the needle through the tunnel.

Suggestion #2

In an alternate technique, attaching the medial and lateral muscles prior to complete implant insertion in the orbit (Fig. 9), the implant tunnels become easier to approach from the side.

Once medial and lateral muscles are attached, the implant should be carefully worked deep into the orbit using a Freer elevator or a small malleable retractor. Adequate deep implant placement is considered important to avoid excessive tension on the overlying Tenon's capsule and conjunctival closures.

After placement is achieved deep into the orbit, the inferior and superior rectus muscles can be sutured to the previously attached medial and lateral muscles (Fig. 10).

Suggestion #3

Placing the suture needle into a needle holder or mosquito forceps so that it points away from the tip of the needle holder (Fig. 11) may also help facilitate placement of the suture, since the direction of the needle will be vertical with respect to the implant and its tunnels, as the holder is tilted to go into the orbit (Fig. 12).

Suggestion #4

Suture placement may also be facilitated by turning the needle holder as much as necessary immediately upon entering a suture hole. This keeps the point of the needle near the top of the tunnel making a shallow pass through it and allowing the needle point to emerge from the exit hole. If the path of the needle is directed too deep in the tunnel, the point will strike the back wall of the tunnel or it will not appear at the exit hole prior to full insertion at tunnel entry. The underlying concept is to "follow the curve of the needle" as the tunnels are designed to accept the curved needle commonly used to attach the extra-ocular muscles.

Techniques for suture placement may be practiced with a non-sterile implant by forming an "eye socket" over the implant using the base of the thumb and forefinger.

Caution:

This document is not intended to be a comprehensive document on surgical techniques. Surgeons should utilize appropriate surgical techniques for which they were trained, and rely upon their clinical experience and judgment, to determine the procedures they will utilize to implant any MEDPOR Surgical Implant. Successful implantations are technique sensitive. Sound surgical judgment should be used in the selection, shaping, handling, and implantation of all MEDPOR shapes.

Additional lanugage translations are available at www.porexsurgical.com.


Fig. 10


Fig. 11


Fig. 12

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
Fig. 11


Fig. 12

Woog, John J., Dresner, Steven, "The Smooth Anterior Surface Suture Tunnel Porous Polyethylene Encapsulation Implant", Poster, presented at the American Society of Ophthalmic Plastic and Reconstructive Surgeons Fall Meeting, New Orleans, LA 2001

Porex Surgical, Inc.

15 Dart Road
Newnan, GA 30265-1017 USA
Email: surgical.info@porex.com
Telephone: 1-678-479-1610 • 1-800-521-7321
Fax: 1-678-423-1435
www.porexsurgical.com

Porex Surgical GmbH

Balanstrasse 55
81541 München, Germany
Email: surgical.info2@porexsurgical.com
Telephone: +49 (0) 89 232415-0
Fax: +49 (0) 89 232415-15
www.porexsurgical.de

© 2010, Porex Surgical Inc. ECO-259-08-B

SST, MCOI, PLUS and EZ are trademarks and MEDPOR is a registered trademark of Porex Surgical, Inc.
The POREX SURGICAL service mark is owned by or used under the license and authority of Porex Corporation